
Synonym	Glycerol, 1,2,3-Propanetriol
----------------	-------------------------------------

State of matter	Liquid
origin	from plants (non-animal origin!)
Density	1.221 - 1.231 kg/l (20°C)
n 20°C/D	1.4490 - 1.4550
Melting point	-10°C
Formula	C ₃ H ₈ O ₃
M	92.10 g/mol
CAS-No.:	56-81-5
HS-No.:	29054500
EC-No.:	200-289-5
Storage:	RT
LGK:	10 - 13
Disposal:	1
WGK:	1

Specification

Assay (titr.)	86.0 - 88.0 %
Aldehyde	max. 0.05 %
Water (K.F.)	approx. 13 %
Chloride	max. 0.0001 %
Fe	max. 0.0005 %
Pb	max. 0.0001 %

Order Information

Prod. No.	Description	Quantity
CE9154	Glycerol 87 %	1000 ml