
State of matter**Liquid**

additional product description

20X concentrated aqueous solution

Storage:

RT

LGK:

10 - 13

WGK:

NA

Specification

DNases/RNases/Proteases

not detectable

pH (20°C; adjusted with HCl)

7.0 ± 0.2

Composition:

tri-Sodium Citrate

88.23 g/L (0.3 M)

Sodium chloride

175.32 g/L (3 M)

Order Information

Prod. No.	Description	Quantity
CE9229	SSC Buffer (20X)	2 l