


Synonym	Polyethyleneglycol-mono-[p-(1,1,3,3-tetra Methylbutyl)-Phenyl]-Ether
State of matter	Liquid
CMC (25°C)	0.3 mM
Classification:	nonionic detergent
Formula	C ₃₄ H ₆₂ O ₁₁
Mass	646.85 g/mol
CAS-No.:	9036-19-5
Storage:	RT
LGK:	10 - 13
Hazard pictogram(s)	
Hazard statement(s)	H302-H318-H411
Precautionary statement(s)	P273-P280-P305+P351+P338-P310
Signal word	Danger
WGK:	2
	®trademark of Union Carbide Company
Specification:	
DNases/RNases/Proteases	not detectable
Heavy metals (as Pb)	max. 0.0005 %
Water	max. 1 %
Chloride	max. 0.005 %
Sulfate	max. 0.005 %

Order Information

Prod. No.	Description	Quantity
CE9240	Triton® X-100	500 ml
CE9241	Triton® X-100	1000 ml