

<b>Synonym</b>	<b>Polyoxyethylene 20 Sorbitane mono-Laurate, Polysorbate 20, Polysorbate 20</b>
State of matter	Liquid
CMC (25°C)	5.9 x 10 <sup>-5</sup> mol/L
Classification:	nonionic detergent
Solubility (25°C)	100 g/L (H <sub>2</sub> O)
Boiling point	> 100°C
Formula	C <sub>58</sub> H <sub>114</sub> O <sub>26</sub>
Mass	1227.72 g/mol
CAS-No.:	9005-64-5
HS-No.:	34021300
NLP-No.:	500-018-3
Storage:	RT
LGK:	10 - 13
Disposal:	28
WGK:	1
	®trademark of Union Carbide Company
Specification:	
DNases/RNases/Proteases	not detectable
Density (d 20°C/4°C)	1.095 - 1.105
Hydroxyl-No.	approx. 100
Saponification No.	40 - 50

## Order Information

<b>Prod. No.</b>	<b>Description</b>	<b>Quantity</b>
CE9242	Tween® 20	1l